Vcheck

Canine CRP 2.0

CANINE REACTIVE PROTEIN

For veterinary use only

INTENDED USE

The Vcheck Canine CRP 2.0 is an *in vitro* diagnostic test kit for the quantitative measurement of C-reactive protein(CRP) concentration in canine serum and plasma. In the case of normal individuals, the concentration of CRP is very low in their serum and plasma but it increases when a pathological phenomenon occurs such as infection, tissue damage or inflammatory disorders.

An increased concentration of CRP is a non-specific reaction to various disease states and can be used to indicate and evaluate immunoinflammatory responses in related diseases, infection and tissue damage.

The Vcheck Canine CRP 2.0 can be used to assist in the diagnosis of a disease state. The test kit is for professional use only.

PRINCIPLE

The Vcheck Canine CRP 2.0 Test kit is a fluorescent immunoassay for the quantitative measurement of canine CRP concentration.

The Vcheck Canine CRP 2.0 Test kit uses specific anti-canine CRP antibodies that will bind to canine CRP. When this specimen is applied to the sample hole of the test device, canine CRP in the specimen migrates along the nitrocellulose membrane and forms complexes with the anti-canine CRP antibodies coated on the membrane. As a result, the density of the test line reflects the concentration of canine CRP in the serum or plasma. The BIONOTE Vcheck Analyzer reads the density of this test line and calculates the canine CRP concentration from the calibration curve data. The control line is a reference line which indicates the test has been performed correctly.

MATERIALS PROVIDED

Reagent	10 Tests/Kit
① Vcheck Canine CRP 2.0 Test device	10
② Assay diluent bottle	10
③ Disposable pipette tips	20
④ Instructions for use	1

MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 5 µl pipette
- 100 µl pipette

STORAGE AND STABILITY

- Store the test kit at 1~30°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

PRECAUTIONS

- The test kit is for canine use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the foil pouch.
- Do not reuse test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers, the components in this kit have been quality control tested as a standard batch unit.
- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.

- Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulations.

COLLECTION AND PREPARATION OF SAMPLE

- Canine serum or plasma should be used with this test.
- [Serum]** Collect the whole blood into a blood collection tube that does not contain anticoagulant. Leave to settle for 30 minutes for blood coagulation and then centrifuge to get a serum supernatant. **[Plasma]** Collect the whole blood into a blood collection tube containing anticoagulant (ONLY heparin). Then centrifuge to get plasma supernatant.
- If serum or plasma samples are not tested immediately, they should be refrigerated at 2~8°C and used within 2 weeks. For longer storage, serum or plasma can be frozen (-20°C or colder). Frozen samples should be brought to room temperature(15~30°C) prior to use.

TEST PROCEDURE

Allow all kit components and sample to reach room temperature (15~30°C) prior to testing.

[Coding]

- Turn on V200 Analyzer and select "Standard Test".


- Remove the test device from the foil pouch. Once the "insert device" is displayed in the screen, insert the test device.


[Dilution of sample & Measurement of CRP value]


- Using a 5 µl pipette, draw 5 µl of sample (serum or plasma) and add the sample into an assay diluent bottle (4 ml).


- Using a 100 µl pipette, mix the sample mixture by pipetting for 5~6 times and draw 100 µl of the sample mixture.

- Add the mixed sample(100 µl) into the sample hole of the test device and press the START to initiate testing.

- The V200 analyzer will display the test result on the screen after 5 minutes.
- Remove the test device.

INTERPRETATION OF THE RESULT

- Read the concentration value of canine CRP appearing on the display of the BIONOTE Vcheck Analyzer.
- If "Low" appears on the display, it means the concentration of canine CRP in the specimen is less than 10 mg/L.
- If "High" appears on the display, it means the concentration of canine CRP in the specimen is greater than 200 mg/L.

REFERENCE RANGE

< 20 mg/L	20~30 mg/L	> 30 mg/L
Normal	Equivocal	Abnormal

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. Solution: Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. Solution: Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. Solution: Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. Solution: Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. Solution: Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23g/dL. Solution: This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely Total hemoglobin	The test is invalid. Solution: Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The calibration is overdue. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	Loading a test device that is not supported by the analyzer. Solution: Check whether the test device is manufactured by BioNote, Inc.
Not Supported Device	Internal error has occurred. Solution: If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
EEE	

PERFORMANCE CHARACTERISTIC

- Analytical specificity
No other canine blood components are known to cross react with CRP in the Vcheck Canine CRP Test kit.
- Interference
Measurement was not influenced by interfering substances such as triglyceride, bilirubin.
- Measuring Range
10~200 mg/L of Canine CRP is measurable in serum and plasma. For >200 mg/L results, samples may be diluted with PBS buffer and re-measured for a numerical result. To calculate the final canine CRP concentration, the measured result must be multiplied by the dilution factor.

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